

# Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/004409

International filing date: 14 February 2005 (14.02.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US  
Number: 60/543,880  
Filing date: 12 February 2004 (12.02.2004)

Date of receipt at the International Bureau: 14 March 2005 (14.03.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland  
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1292555

# THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

*March 04, 2005*

**THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.**

**APPLICATION NUMBER: 60/543,880**

**FILING DATE: February 12, 2004**

**RELATED PCT APPLICATION NUMBER: PCT/US05/04409**



Certified by

Under Secretary of Commerce  
for Intellectual Property  
and Director of the United States  
Patent and Trademark Office

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

**PROVISIONAL APPLICATION FOR PATENT COVER SHEET**

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. EE 742523104 US

**INVENTOR(S)**

| Given Name (first and middle [if any]) | Family Name or Surname | Residence<br>(City and either State or Foreign Country) |
|--|------------------------|---|
| George                                 | Tzertzinis             | Cambridge, MA   |
| George                                 | Feehery                | West Newbury, MA  |
| Larry                                  | McReynolds             | Beverly, MA   |

Additional inventors are being named on the 1 separately numbered sheets attached hereto**TITLE OF THE INVENTION (500 characters max)**

Active Heterogeneous siRNA Mixtures

Direct all correspondence to:

**CORRESPONDENCE ADDRESS**
☒ Customer Number: 28986  
 OR

|   |                           |           |              |     |              |
|---|---------------------------|-----------|--------------|-----|--------------|
| <input checked="" type="checkbox"/> Firm or Individual Name | New England Biolabs, Inc. |           |              |     |              |
| Address   | 32 Tozer Road             |           |              |     |              |
| Address   |                           |           |              |     |              |
| City  | Beverly                   | State     | MA           | ZIP | 01915        |
| Country   |                           | Telephone | 978-927-5054 | Fax | 978-927-1705 |

**ENCLOSED APPLICATION PARTS (check all that apply)**

|   |  |
|---|--|
| <input checked="" type="checkbox"/> Specification Number of Pages <u>14</u> | <input type="checkbox"/> CD(s), Number _____   |
| <input type="checkbox"/> Drawing(s) Number of Sheets _____                  | <input type="checkbox"/> Other (specify) _____ |
| <input type="checkbox"/> Application Data Sheet. See 37 CFR 1.76            |  |

**METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT**

|   |   |
|---|---|
| <input checked="" type="checkbox"/> Applicant claims small entity status. See 37 CFR 1.27.  | FILING FEE<br>Amount (\$)<br><br><b>80.00</b> |
| <input checked="" type="checkbox"/> A check or money order is enclosed to cover the filing fees.  |   |
| <input type="checkbox"/> The Director is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number: _____ |   |
| <input type="checkbox"/> Payment by credit card. Form PTO-2038 is attached.   |   |

The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.

|  |
|--|
| <input checked="" type="checkbox"/> No.  |
| <input type="checkbox"/> Yes, the name of the U.S. Government agency and the Government contract number are: _____ |

[Page 1 of 2]

Respectfully submitted,

SIGNATURE

TYPED or PRINTED NAME Harriet M. StrimpelTELEPHONE 978-927-5054Date 02/12/04REGISTRATION NO. 37,008

(If appropriate)

Docket Number: NEB-239**USE ONLY FOR FILING A PROVISIONAL APPLICATION FOR PATENT**

This collection of information is required by 37 CFR 1.51. The information is required to obtain or retain a benefit by the public which is to file (and by the USPTO to process) an application. Confidentiality is governed by 35 U.S.C. 122 and 37 CFR 1.14. This collection is estimated to take 8 hours to complete, including gathering, preparing, and submitting the completed application form to the USPTO. Time will vary depending upon the individual case. Any comments on the amount of time you require to complete this form and/or suggestions for reducing this burden should be sent to the Chief Information Officer, U.S. Patent and Trademark Office, U.S. Department of Commerce, P.O. Box 1450, Alexandria, VA 22313-1450. DO NOT SEND FEES OR COMPLETED FORMS TO THIS ADDRESS. SEND TO: Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

If you need assistance in completing the form, call 1-800-PTO-9199 and select option 2.

021204

15535 U.S. PTO  
60/543880

021204

**PROVISIONAL APPLICATION COVER SHEET**  
**Additional Page**

PTO/SB/16 (08-03)

Approved for use through 07/31/2006. OMB 0651-0032

U.S. Patent and Trademark Office; U.S. DEPARTMENT OF COMMERCE

Under the Paperwork Reduction Act of 1995, no persons are required to respond to a collection of information unless it displays a valid OMB control number.

Docket Number

NEB-239

**INVENTOR(S)/APPLICANT(S)**

| Given Name (first and middle [if any]) | Family or Surname | Residence<br>(City and either State or Foreign Country) |
|--|-------------------|---|
| Derek                                  | Robinson          | Boxford, MA   |
| Sriharsa                               | Pradhan           | Wenham, MA  |

[Page 2 of 2]

Number 1 of 1

**WARNING:** Information on this form may become public. Credit card information should not be included on this form. Provide credit card information and authorization on PTO-2038.

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re application of: Tzertzinis et al.

Application No.: not yet assigned  
Filed: herewith  
For: Active Heterogeneous siRNA Mixtures

Group No.: N/A  
Examiner: N/A

Mail Stop Provisional Application  
Commissioner for Patents  
P.O. Box 1450  
Alexandria, VA 22313-1450

EXPRESS MAIL CERTIFICATE


"Express Mail" label number EE 742523104 US  
Date of Deposit: February 12, 2004

I hereby state that the following *attached* paper or fee

Provisional patent application and related papers

is being deposited with the United States Postal Service "Express Mail Post Office to Addressee" service under 37 C.F.R. § 1.10, on the date indicated above and is addressed to the Mail Stop Provisional Application, Commissioner for Patents, P.O. Box 1450, Alexandria, VA 22313-1450.

Leslie Goldberg



Signature of person mailing paper or fee

Docket No. : NEB-239

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE  
APPLICATION FOR UNITED STATES LETTERS PATENT

Inventors: George Tzertzinis  
George Feehery  
Larry McReynolds  
Derek Robinson  
Sriharsa Pradhan

Title: ACTIVE HETEROGENEOUS siRNA MIXTURES

Attorney: Harriet M. Strimpel, D. Phil.  
Patent Counsel  
New England Biolabs, Inc.  
32 Tozer Road  
Beverly, MA 01915

Customer No. 28986

## **ACTIVE HETEROGENEOUS siRNA MIXTURES**

**G.Tzertzinis, G. Feehery, L. McReynolds, D. Robinson, S.  
Pradhan**

5

### **THIS IS A PROVISIONAL APPLICATION**

#### **BACKGROUND OF THE INVENTION**

There is an increased demand for the generation of gene-specific  
10 inactivation reagents for functional genomics, for genes that are  
potential drug targets, and biological pathway investigations. Gene  
silencing double-stranded RNA fragments (siRNAs) have been used in  
many cases to specifically inactivate selected targets in mammalian  
cells. Major obstacles in such endeavors are the costly design  
15 methodology and discovery of effective siRNAs since individual siRNAs  
vary greatly in their efficiency for gene silencing due to a phenomenon  
generally termed position dependence effect. This requires the  
synthesis and testing of multiple siRNAs to identify the most effective  
ones, a process known as siRNA validation. These stages in the  
20 production of effective siRNAs ready for use raises several fold the cost  
and time required for obtaining a specific reagent for each target.  
Additionally chemical synthesis of RNA is inherently expensive.

Standard methodologies require the use of 20-150 nM of siRNA in  
25 each transfection for appreciable effects. Recent reports have,  
however, demonstrated that concentrations higher than 25 nM  
transfected into mammalian cells lead to non-specific effects known as  
“off-target” effects (Semizarov et al.) which illustrates the necessity of  
validating each siRNA.

This application relates to the generation of gene-specific mixtures of siRNA that can be used as ready-for-use reagents for gene inactivation studies.

## 5                    **Description of Embodiments**

1. Cloned DNA suitable for generating double-stranded RNA substrates for heterogeneous siRNA mixtures

10        Segments of mRNA sequence were selected for an identified target preferably after a sequence comparison within a gene database. An algorithm was developed to scan for candidate gene sequences using a window of 16-21 bp and providing matching sequences in the database. (for example, the UNIGENE database). Regions of the target  
15        sequence showing absence of hits to other targets were preferably selected.

          Although there is no particular limitation on segment size, in the examples herein segments were selected in the range of about 150-  
20        1000 bp for example, 200-400 bp long.

          The selected segment was amplified with PCR primers made according to standard protocols for PCR primer design. For example, primers may include a T7 promoter sequence at the 5' end. Other  
25        sequences can be used in place of the T7 promoter to facilitate cloning to one of the double T7 promoter vectors (Litmus 28i, Litmus 38i, Litmus-U from NEB).



For Litmus U, the following primer sequences were used: gggaaagu and ggagacau, where u stands for uracil. After the PCR reaction, the amplified DNA product was cloned directly in Litmus U using the USER protocol (NEB). The cloned fragments were used for the production of dsRNA using HiScribe (NEB). dsRNA was prepared from DNA (see above for example) using HiScribe (NEB). The reaction mixture was incubated at 42 °C for three hours and the dsRNA was phenol extracted, and dialyzed.

## 2. Generation and purification of hsiRNA

21-22 bp hsiRNA was typically produced from an RNase III digestion in the presence of manganese ions (see US patent application 10/622,240 incorporated by reference) and (Shortcut RNA kit from NEB). ShortCut siRNA is a highly potent mixture of 21-22 bp dsRNA that is processed from a large gene-specific dsRNA (150-1500 bp) by RNase III in the presence of manganese buffer, the sequences of which collectively span the entire target RNA. The large dsRNA construct has been designed to be uniquely representative of the target gene with as little homology as possible to other areas of the genome.

The material obtained from the cleavage reaction was purified on a Pharmacia Source 15Q packed in an HR16/10 in a salt gradient using the following buffers:

Buffer A: 0 M NaCl, 20 mM Tris-Hcl pH 7.5 (25° C), 0.5 mM EDTA.

Buffer B: 1 M NaCl, 20 mM Tris-Hcl pH 7.5 (25° C), 0.5 mM EDTA.

The column was run on a Pharmacia AKTA FPLC system using the following program parameters:

- 5    Flow rate: 6 ml/min.  
     Start Concentration B: 0%.  
     Equilibration: 10 column volumes.  
     Load volume: ~100 ml.  
     Wash: 2 column volumes buffer A.
- 10   Fraction Size: 6 ml.  
     Target Concentration: 100%B.  
     Length of gradient: 20 column volumes.  
     Detection wavelength: 260nm
- 15   Pooling strategy/Source 15Q resolution

- The Source Q resolves the major ds RNA digestion product (~21-22 base pair long) from smaller RNA fragments and RNase III which elute from the column at lower salt prior to the collected peak, and from
- 20   large dsRNA and substrate DNA which elute after the major peak at higher salt, observed by monitoring the flow conductivity. The siRNA containing fractions were pooled so as to avoid any larger size material, and any smaller RNA, while capturing the majority of the 21 bp siRNA. Fractions analyzed by non-denaturing PAGE confirmed this
  - 25   separation profile.

Examples of siRNA mixture preparations purified on the Source 15Q showed that the fragments of different hsiRNA were eluted at similar

salt concentrations making it possible to standardize the protocol for any desired fragment.

- 5 a. GFP siRNA fractions 30,31 eluted at 47.6% B and a conductivity of 43%.
- b. Luciferase siRNA fractions 30-32 eluted at 48.3% B and a conductivity of 43.6%.
- c. Litmus 28i polylinker siRNA fractions 29,30 eluted at 46%B and a conductivity of 41.4%.
- 10 d. Erk2 siRNA fractions 29-32 eluted at 47.6%B and a conductivity of 41.5%.

The Source 15Q pooled fractions of siRNA mixture from 3-4 fractions (usually 18-24 mls) were then dialyzed overnight against 2 liters of  
15 Storage Buffer ( 20 mM KCl, 10mM HEPES pH 7.0 (at 25<sup>0</sup> C), 0.5 mM EDTA, made with milli-Q or equivalent water). The dialyzed siRNA was frozen at -20 °C. Sterilize using passage through 0.2 µm filter and adjust the concentration to 150 ng/µl by dilution with sterile storage buffer before freezing for storage.

20

### 3. Endotoxin removal from heterogeneous siRNA mixture

For some applications it may be desirable to complete deplete the siRNA mixtures from any endotoxin (LPS) carried over from previous  
25 manipulations.

Endotoxin was measured using an LAL pyrochrome kit from Associates of Cape Cod Inc, Falmouth MA CAT# C0180. The endpoint method

listed in the product literature was used. The endotoxin was measured in Endotoxin Units (EU).

5 Examples of endotoxin levels measured in Source 15Q pools before purification are provided below:

- a. AKT siRNA at approximately 0.5 mg/mL gave 16 EU/ml.
- b. Luciferase siRNA gave 10.4 EU/ml.

10 Endotoxin levels can be significantly reduced by using a Pharmacia Source RPC directly after the Source 15Q protocol (3) to levels below 1EU.

The Source 15Q pool was loaded on a 3 ml Source RPC column.

15

Buffer A: 35 mM triethylamine (pH 7.0 with acetic acid at 25<sup>0</sup> C), 2% acetonitrile made with Milli-Q or equivalent water.

Buffer B: 100% acetonitrile.

20 The column was run on a Pharmacia AKTA FPLC system using the following program parameters:

Flow Rate: 2 ml/min.

Start Concentration B: 0%.

Equilibration: 5 column volumes A.

25 Load volume: ~50-100 ml.

Wash: 2 column volumes A.

Fraction Size: 1 ml.

Target Concentration: 20% B.

Length of gradient: 15 column volumes.

This protocol removes any contaminating endotoxin from the siRNA mixture which elutes at approximately 50% buffer B.

Examples of siRNA preparations purified on the Source RPC.

5

- a. AKT siRNA fractions 27,28 eluted at 10.7% B.
- b. Luciferase siRNA fractions 24-27 eluted at 10.6% B.

Fractions were transferred to 1.5 ml micro-centrifuge tubes and dried overnight in a spin vacuum without heat. Pellets were hydrated with storage buffer at room temperature.

10

Examples of endotoxin levels after the RPC purification:

- 15 a. AKT siRNA fraction 27 contained 0.18 EU/ml (800X reduction),  
fraction 28 contained 0.72 EU/ml (200X reduction).
- b. Luciferase siRNA pool contained 0.062 EU/ml (168X reduction).

20

Both the Luciferase source Q- purified and RPC-purified siRNA mixtures are able to knock down by over 90% Luciferase expression transfected with siRNA mixtures at 1 nM in COS cells.

#### 4. Preparation of a Kit

Applications of the kit include

25

- Gene silencing
- Target validation

The hsiRNA was formulated in: 20 mM KCl, 10 mM Na-HEPES (pH 7.0), 0.5 mM EDTA and was free of contaminating large molecular weight dsRNA, ssRNA, DNA, and protein.

The kit further optionally contains a transfection mixture (Transit-IT-  
5 TKO (Mirus Corp))

## 5. Transfection Protocol for hsiRNA into cells:

- A very small amount of the hsiRNA is sufficient for effective  
10 silencing as compared to single sequence siRNA's. A starting concentration of 20 nM can be used which corresponds to 1  $\mu$ l (10 pmol, 150 ng) of hsiRNA in 0.5 ml of transfection media.
- 15 a. Plate cells the day before so that the cell density is 40-60% confluent at the time of transfection.
  - b. Mix an appropriate amount of hsiRNA transfection reagent with serum free medium. Incubate at room temperature for 10-20 minutes.
  - 20 c. Add an appropriate volume of hsiRNA mix (see table) to the diluted transfection reagent and incubate 10-20 minutes at room temperature to form the transfection complexes.
  - d. Dilute the complex with complete medium to the desired final culture volume for the plate size used (see table).
  - 25 e. Aspirate the medium from the cell plate and replace with the diluted transfection complex.

Incubate cells 24-48 hours before analysis.

Examples for one transfection per well of the indicated size plates are shown in the table below. In the volumes shown the final siRNA mixture concentration is 20 nM.

| Plate size           | 6               | 12             | 24          | 96             |
|----------------------|-----------------|----------------|-------------|----------------|
| transfection reagent | 6-10<br>$\mu$ l | 4-6<br>$\mu$ l | 2-4 $\mu$ l | 1-2 $\mu$ l    |
| Serum free medium    | 200<br>$\mu$ l  | 100<br>$\mu$ l | 50 $\mu$ l  | 25 $\mu$ l     |
| siRNA mix            | 4<br>$\mu$ l    | 2<br>$\mu$ l   | 1 $\mu$ l   | 0.1<br>$\mu$ l |
| Complete medium      | 800<br>$\mu$ l  | 650<br>$\mu$ l | 450 $\mu$ l | 75 $\mu$ l     |
| Final volume         | 1000<br>$\mu$ l | 750<br>$\mu$ l | 500 $\mu$ l | 100 $\mu$ l    |

5

## 6. Examples of Active Mixtures

10


Name: CREB ShortCut siRNA mix  
 Concentration: 10  $\mu$ M (150 ng/ $\mu$ l)  
 Size: 15  $\mu$ g

15

Description: A heterogeneous mixture of 21-22 bp short interfering RNAs (siRNA) that induces effective silencing (RNAi) of the endogenous transcription factor CREB at concentrations of 20 nM and below in mammalian cell lines. CREB is a member of the leucine zipper

family of DNA binding proteins, binds as a homodimer to the cAMP-responsive element (CRE) and activates transcription in response to a variety of extracellular signals.

## 5 Western Blot:



- 15 Western blot analysis of extracts from HeLa cells transfected with CREB ShortCut siRNA mix (+) or control siRNA (-). An antibody to the targeted protein CREB confirms silencing of protein expression, while a antibody against non-targeted AKT is used to control protein loading and to confirm siRNA specificity.

20

Source: A 360 bp DNA template derived from a CREB Mouse cDNA construct (coordinates 247-607, accession number m34356) was transcribed in vitro by T7 RNA polymerase using HiScribe to create double-stranded RNA (dsRNA). The dsRNA was processed by RNase III in the presence of manganese buffer (Shortcut) to produce a mixture of 21-22 bp siRNAs, and purified by column chromatography.

25

Name: p38 MAPK1 hsiRNA



Concentration: 10  $\mu$ M (150 ng/ $\mu$ l)

Size: 15  $\mu$ g


Description: A heterogeneous mixture of 21-22 bp short interfering  
5 RNAs (siRNA) that induces effective silencing (RNAi) of the MAP kinase  
P38 at concentrations of 20 nM or less in mammalian cell lines.

Supplied in: 20 mM KCl, 10 mM Na-HEPES (pH 7.0), 0.5 mM EDTA.

Source: A 418 bp DNA template derived from a Human p38 cDNA  
construct (coordinates 10-419, accession number L35253) is  
10 transcribed by T7 RNA polymerase to create double-stranded RNA  
(dsRNA). RNase III cleaves the dsRNA in the presence of manganese  
buffer to 21-22 bp siRNA.

Western Blot:

15



20

p38 (target)



ATP citrate lyase (control)

25

- + ShortCut p38 siRNA (6.6 nM)

Western blot analysis of extracts from HeLa cells transfected with p38 ShortCut siRNA mix (+) or non-targeted (-) control siRNA. An antibody to the targeted protein p38 confirms silencing of protein expression, while a non-targeted antibody against cs is used to control protein loading and to confirm siRNA specificity.

Name: p42/44 MAPK1 (ERK2) ShortCut siRNA mix


Concentration: 10  $\mu$ M (150 ng/ $\mu$ l)

Size: 15  $\mu$ g

Description: A heterogeneous mixture of 21-22 bp short interfering RNAs (siRNA) that induces effective silencing (RNAi) of the endogenous transcription factor ERK2 at concentrations of 20 nM and below in mammalian cell lines.

Source: A 283 bp DNA template derived from a Human Erk1 cDNA construct (coordinates 667-950, accession number NM\_002745) is transcribed by T7 RNA polymerase to create double-stranded RNA (dsRNA). RNase III cleaves the dsRNA in the presence of manganese buffer to 21-22 bp siRNA.

Western Blot:



Western blot analysis of extracts from HeLa cells transfected with  
15 ERK2 ShortCut siRNA mix (+) or non-targeted (-) control siRNA. An  
antibody to the targeted protein ERK2 confirms silencing of protein  
expression, while a non-targeted antibody against Dynamin is used to  
control protein loading and to confirm siRNA specificity.

Table 1

| CATEGORY       | TARGET            | <b><i>Acc. number</i></b> | Coordinates |
|----------------|-------------------|---------------------------|-------------|
| kinases        | Akt1              | NM_005163                 | 199-657     |
|                | Erk2              | NM_002745                 | 660-940     |
|                | MSK               | AF074393                  | 282-736     |
|                | p38               | L35253                    | 10-419      |
|                | IRS1              | NM005544                  | 1026-1713   |
|                | PKR               | M35663                    | 999-1499    |
|                | PTEN              | NM_000314                 | 1019-1445   |
| transcription  | CREB              | M34356                    | 247-601     |
| Nuc. signaling | ERa               | NM_000125                 | 369-905     |
|                | ERb               | NM_001437                 | 587-1240    |
|                | DAX               | NM_000475                 | 1-249       |
|                | p53               | NM_000546                 | 717-915     |
|                | DNMT1             | X69632-G-BPR2             | 2124-3235   |
|                | DnMT3B            | AF331857                  | 1150-1545   |
|                | DnMT3A            | X63692.gb-pr2             | 1547-2388   |
|                | TRIP              | L38810                    | 1-445       |
|                | Rb                | m15400.gb_pr1             | 2239-2755   |
|                | MeCP2             | af030876.gb_pr            | 699-1011    |
| Other          | caspase3          | p42574                    | 1063-1496   |
|                | La                | NM_003142                 | 316-631     |
|                | FURIN             | NM002569                  | 1781-1990   |
| Controls,      | Lit28i polylinker | NEB#N3528S                | 2465-2600   |
| gen. use       | EGFP              | U55763                    | 596-1322    |
|                | RFP               | AF272711                  | 152-632     |
|                | FfLUC             | U47295                    | 747-1757    |
|                | Renilla           | AF264722                  | 3673-3951   |

Table 1 lists a series of target genes for which hsiRNA fragments have  
5 been or are being prepared from dsRNA having a sequence  
corresponding to the coordinates for the gene (cDNA) sequence  
contained in accession number of GenBank given above.